

Agarose LE

Cat. No. 32034 / 500g

The first and finest name in agarose

- ❖ Analytical grade(Routine analysis)
- ❖ Gel strength(1%) > 1,200g/cm²
- ❖ GellingTemp range (1.5%) 36 °C ± 1.5 °C
- ❖ MeltingTemp ≥ 90 °C
- ❖ DNase/RNase Non detected

Applications

Analytical electrophoresis of DNA and RNA ≥ 100 bp
Northern / Southern blotting of DNA and RNA

