

MycoClean™ Mycoplasma Prevention Spray

Cat. No. 21083 / 200ml

MycoClean™ Mycoplasma Prevention Spray is an effective solution for reliable disinfection of laboratory surfaces and apparatus including clean benches, incubators, work benches, cell storage boxes and liquid nitrogen containers. MycoClean™ shows an excellent effectiveness against mycoplasma. Contamination-free and clean surfaces can be achieved in a single application.

- ❖ Eradicates mycoplasma
- ❖ Non-corrosive and non-carcinogenic
- ❖ Convenient and stable
- ❖ Ready-to-use

FEATURES

MycoClean® originated from a natural substance. eradicates mycoplasma by forming pore to the cell membrane, using osmotic pressure.

INSTRUCTION FOR USE

MycoClean™ Mycoplasma Prevention Spray is a ready-to-use solution supplied in a refillable spray bottle. The use of MycoClean™ is fast, easy and effective to maintain a clean work area.

STORAGE CONDITION

Store at room temperature. and then stable for more than 2 years.

INSTRUCTION FOR USE

For *in vitro* research use on

